

**Product Name :**  
Classic Eye Model

**Product Code :**  
ANATOMYSUK-AS06LAB004D



**Description :**

Classic Eye Model

**Technical Specification :**

Classic human eye model is 3X life-size and splits into 6 parts.

Removable parts include: both halves of the sclera, both halves of the choroid, eye lens, and vitreous humour.

Mounted on a circular base.

Great model to teach and learn the anatomy of the eye.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265