

Product Name :
Theodolite Model

Product Code :
MATHSUK-AS05LAB002A



Description :

Theodolite Model

Technical Specification :

A theodolite is an instrument for measuring both horizontal and vertical angles, as used in different types of works as triangulation, prolonging, computation of elevation and depression of distant and near.

It consists of movable telescope mounted on the horizontal and vertical axes. Both the axes of theodolite are equipped with graduated circles.

A rotating telescope for measuring horizontal and vertical angle which is conceptually it is a great way to understand the mechanics of a theodolite.

Inculcates a spirit of experimenting among the students.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265