

**Product Name :**  
Coaxial Binocular Microscope

**Product Code :**  
MICROSCOPESUK-AS11LAB005A



**Description :**

Coaxial Binocular Microscope

**Technical Specification :**

High quality Binocular Microscope for School, Colleges and Research Institutes.

Provided with a Binocular observation head, inclined at 45 degree, with interpupillary distance 55mm to 75 mm & diopter adjustments.

Rectangular 115 x 125 mm with co-axial mechanical stage.

Box type base with 6v-20w, variable halogen illumination.

Sub stage movable by rack and pinion Abbe Condenser of 1.25 N.A. with iris diaphragm and filter

Achromatic objectives:4x, 10x, 40x, 100x(SL).

Huygenian eyepieces:5x (pair).

Eyepieces:10x(pair) WF.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265