

Product Name :
Fisher Clamp (Routine Improved Model)

Product Code :
CHEMISTRYSUK-AS03LAB002A



Description :

Fisher Clamp (Routine Improved Model)

Technical Specification :

Double burette clamp with all aluminium casted mounts on rod of 1/4" to 1/2" in diameter, two burettes can be securely held in a vertical position whereas still permitting full read of the burette scale.

Clamp includes a sleek enamel end with rubber finger grippers.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com