

Product Name :
Fisher Clamp (New Improved Model)

Product Code :
CHEMISTRYSUK-AS03LAB003A



Description :

Fisher Clamp (New Improved Model)

Technical Specification :

Fisher sort clamps, for holding burettes, complete body of formed metal creating it light-weight, sturdy and durable.

The special style provides vice-like grip to the measuring system for a positive clamping.

Given a spring loaded positive grip and unharness.

Bosshead clamp secures it firmly at desired position on the rod of regarding 13mm diameter.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265