

**Product Name :**  
Adjustable Clamp Holder

**Product Code :**  
CHEMISTRYSUK-AS03LAB004A



**Description :**

Adjustable Clamp Holder

**Technical Specification :**

For rods up to 5/8" in diameter.

Clamp holder that is adjustable to swivel and lock to any angle.

Enamel finish, tightened with a plastic thumb screws.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

