


Product Name :
Balance Triple Beam

Product Code :
CHEMISTRYSUUK-AS03LAB001C



Description :

Balance Triple Beam

Technical Specification :

A single pan low balance with three tiered beams which can weigh upto 2610g when used with supplementary masses, the balances has three notched weighing beam with center indicating sliding masses giving total capacity of 610g.

Additional masses are hung on the end of the beam increasing the capacity upto 2610g.

And they can be stored in the recessed base when not required

Magnetic Damping.

Beams: 0 to 500g x 100g

0 to 100g x 10g

0 to 10g x 0.1g

Supplementary masses: 2 x 1000g

Pan, stainless steel dia.: 150 mm.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265