

**Product Name :**  
Carboy With Stop Cock

**Product Code :**  
CHEMISTRYSUK-AS03LAB004D



**Description :**

Carboy With Stop Cock

**Technical Specification :**

These carboys are an excellent replacement for glass carboys and are used for storing and dispensing distilled water & other laboratory solutions.

These polypropylene Carboys are autoclavable, dependable and durable.

These are also equipped with a leak proof spigot screwed to the body.

**Material:** Polypropylene

**Capacity:** 10 Ltrs., 20 Ltrs.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265