

Product Name :
Aspirator Bottles

Product Code :
CHEMISTRYSUK-AS03LAB005D



Description :

Aspirator Bottles

Technical Specification :

Aspirator Bottles are made of Polypropylene and are therefore much lighter than glass. These are autoclavable, durable and dependable. The leakproof spigot screwed to the bottle body requires only a half turn to open or close. These bottles are used for storing or dispensing distilled water & other laboratory solutions. Threads should be totally disengaged before autoclaving.

Material: Polypropylene

Capacity: 5 Ltrs., 10 Ltrs., 20 Ltrs.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265