

**Product Name :**  
Bottle Reagent, W.M. Glass

**Product Code :**  
CHEMISTRYSUK-AS03LAB001E



**Description :**

Bottle Reagent, W.M. Glass

**Technical Specification :**

Wide mouth with dustproof flat stopper of Glass/Polythene.

**Capacity:-**

60ml

125ml

250ml

500ml

1000ml



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265