

Product Name :
Bottle Reagent, N.M. Glass Amber Colour

Product Code :
CHEMISTRYSUK-AS03LAB002E



Description :

Bottle Reagent, N.M. Glass Amber Colour

Technical Specification :

Narrow mouth with dustproof flat stopper of Glass/Polythene.

Capacity:-

30ml

60ml

125ml

250ml

500ml

1000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265