

Product Name :
Bottle Reagent, W.M. Glass Amber Colour

Product Code :
CHEMISTRYSUK-AS03LAB003E



Description :

Bottle Reagent, W.M. Glass Amber Colour

Technical Specification :

Soda glass wide mouth with dustproof flat stopper

Capacity:-

30ml

60ml

125ml

250ml

500ml

1000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265