

**Product Name :**  
Polypropylene Beaker

**Product Code :**  
CHEMISTRYSUK-AS03LAB001F



**Description :**

Polypropylene Beaker

**Technical Specification :**

Beaker, moulded in Polypropylene, have excellent clarity and very good chemical resistance.

Pouring a solution becomes very convenient because of its tapered spout.

These autoclavable Beakers have prominently raised graduations on the outer surface.

**Capacity:-**

50ml

100ml

250ml

500ml

1000ml

2000ml



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265