


Product Name :
Measuring Jugs

Product Code :
CHEMISTRYUK-AS03LAB005F



Description :

Measuring Jugs

Technical Specification :

These handy jugs, moulded in Polypropylene, are clear, autoclavable and have good chemical resistance.

These jugs are provided with a handle for easy and convenient lifting and have raised graduations for easy reading.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com