

Product Name :
Wide Mouth Reagent Bottles

Product Code :
CHEMISTRYSUK-AS03LAB002G



Description :

Wide Mouth Reagent Bottles

Technical Specification :

These Wide Mouth Reagent Bottles are available in Polypropylene. The PP bottles are rigid, translucent and have excellent contact clarity. These bottles can be autoclaved filled as well as empty. The Polypropylene Cap has a built in seal ring to make these bottles leakproof.

Material: Polypropylene

Capacity:-

60ml

125ml

250ml

500ml

1000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265