

**Product Name :**  
Burette

**Product Code :**  
CHEMISTRYUK-AS03LAB002H



**Description :**

Burette

**Technical Specification :**

This Burette has a transparent acrylic body and a leakproof stopcock made fitted with a self lubricating PTFE plug.

Each Burette is individually calibrated to give maximum accuracy.

Acrylic body is unaffected by dilute mineral acids and bases except Hydrofloric Acid & Ammonium Hydroxide.

The pH range of 1-14 is ideal and can be used up to 5% wt/wt of acids or 30% bases.

**Material:** Polymethyl Methacrylate/TPX

**Capacity:-**

25ml

50ml

100ml



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265