

**Product Name :**  
Measuring Cylinder Hexagonal

**Product Code :**  
CHEMISTRYUK-AS03LAB003J



**Description :**

Measuring Cylinder Hexagonal

**Technical Specification :**

Measuring Cylinders are now available with hexagonal base also.

The hexagonal base provides excellent stability to the cylinder where as the tapered spout works really well while pouring any reagent.

**Material:** Polypropylene

**Capacity:-**

10ml

25ml

50ml

100ml

250ml

500ml

1000ml



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265