

Product Name :
Borosilicate Glass Filtering Funnel

Product Code :
CHEMISTRYUK-AS03LAB005L



Description :

Borosilicate Glass Filtering Funnel

Technical Specification :

Made of short stem, 60° conical shape.

Capacity:-

50ml

62ml

75ml

100ml

150ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265