


Product Name :
Retort Stand

Product Code :
CHEMISTRYSUK-AS03LAB001M



Description :

Retort Stand

Technical Specification :

Cast iron rectangular base with stove hammertone painted with hole.
12.5mm dia M.S. Rod chrome plated one end half round and other end
machine threaded with washer and hexagonal nut.

Base of
Retort
Stand

(mmHg)		
1500	1500 Kg	
2000	2000 Kg	
2500	2500 Kg	
2750	2750 Kg	
3000	3000 Kg	

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific