

**Product Name :**  
Test Tube Stand 3 Tier

**Product Code :**  
CHEMISTRYUK-AS03LAB003M



**Description :**

Test Tube Stand 3 Tier

**Technical Specification :**

The PC rack has transparent upper & middle plates where as in PP racks, every part is colored & opaque. These ingenious racks with a 3-tier design provide clear view of the tube contents & convenient insertion & removal of tubes. These auto clavable Racks can withstand sub-freezing temperatures.

Material : Polypropylene.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

