

Product Name :
Funnel Holders

Product Code :
CHEMISTRYSUK-AS03LAB005M



Description :

Funnel Holders

Technical Specification :

A Non-corrosive funnel holders have outer ring to hold funnels with 3"-6" dia while the removable inner lid holds smaller funnels with 1"-3" dia.

Available :

- Single Funnel Holder.
- Double Funnel Holder.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265