

Product Name :
Bottle Reagent, N.M. Glass Amber Color

Product Code :
CHEMISTRYSUK-AS03LAB001N



Description :

Bottle Reagent, N.M. Glass Amber Color

Technical Specification :

With interchangeable stopper.

Material : Borosilicate glass.

Heat Resistant.

Narrow mouth.

Chemical Resistant.

Amber Color.

Capacity:

30ml

60ml

125ml

250ml

500ml

1000ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265