

Product Name :
Bottle Reagent, W.M. Glass Amber Color

Product Code :
CHEMISTRYSUK-AS03LAB002N



Description :

Bottle Reagent, W.M. Glass Amber Color

Technical Specification :

Soda glass.

Material : Borosilicate glass.

Heat Resistant.

Chemical Resistant.

Amber Color.

Wide mouth with dust proof flat stopper.

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific