


**Product Name :**  
Fisher Clamp, Metal Die Cast

**Product Code :**  
CHEMISTRYSUK-AS03LAB004P



**Description :**

Fisher Clamp, Metal Die Cast

**Technical Specification :**

The special design provides vice like grip to the burette for positive clamping.

Bosshead clamp secures it firmly at desired position on the rod of about 13mm diameter.

Fisher type clamps, for holding burettes, complete body of die cast metal making it light weight, durable and sturdy.

Provided with a spring loaded positive grip and release.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

