

**Product Name :**  
Oak Ridge Centrifuge Tubes

**Product Code :**  
CHEMISTRYSUK-AS03LAB005P



**Description :**

Oak Ridge Centrifuge Tubes

**Technical Specification :**

Oak Ridge Centrifuge Tubes, made of Polypropylene, have excellent contact clarity and can be used upto 50000 x G.

**Capacity:** 30ml, 50ml.

The screw cap should not be used below 4°C.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265