

Product Name :
MCT Box

Product Code :
CHEMISTRYUK-AS03LAB001Q



Description :

MCT Box

Technical Specification :

MCT Boxes are moulded in polycarbonate and are therefore strong as well as autoclavable.

The transparent cover provides a good view of the contents inside and is numbered for better inventory control.

One size can hold 64 MCTs of 1.5ml and the other can hold 100 MCTs of 0.5ml.

These boxes are designed to work at temperatures ranging from -190°C to + 121°C.

Capacity:-

Rack for 64 MCTs of 1.5ml

Rack for 100 MCTs of 0.5ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265