

**Product Name :**  
Sample Container (Screw Type)

**Product Code :**  
CHEMISTRYUK-AS03LAB002Q



**Description :**

Sample Container (Screw Type)

**Technical Specification :**

These thin walled, disposable containers are easily screw capped with less than 1/4th of a rotation.

This is perhaps the most ideal option available for carrying samples from one place to another safe and sound.

Sample Containers are made of polypropylene.

As one screw caps the container, the top of it sets in a domed ring in the cap making it completely leakproof.

**Capacity:** 25ml, 30ml

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**