

**Product Name :**  
Urine Container

**Product Code :**  
CHEMISTRYSUK-AS03LAB003Q



**Description :**

Urine Container

**Technical Specification :**

Moulded in polypropylene these container have excellent contact clarity and the cap has built-in seal that makes it completely leakproof.

These containers can be used for collecting a wide range of samples of sputum, urine etc.

**Capacity:** 30ml, 60ml.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265