


Product Name :
Test Tube With Screw Cap

Product Code :
CHEMISTRYSUUK-AS03LAB006Q



Description :

Test Tube With Screw Cap

Technical Specification :

The cap is screwed on to the tube & hence this tube can be filled up to the brim.

The polypropylene tubes provide excellent option for Ria, coagulation and bacteriology etc.

Capacity: 16 x 100mm

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

