


**Product Name :**  
Test Tube Cap

**Product Code :**  
CHEMISTRYSUK-AS03LAB001R



**Description :**

Test Tube Cap

**Technical Specification :**

The grooves inside the cap tighten the fitting on the tube while the outer grooves would hold identification labels.

This transparent polypropylene cap for test tubes proves to be an excellent protection to the tubes.

The built-in seal ring would set in the tubes making the combination a leakproof one.

**Capacity:** Cap for Test Tube Dia 25mm



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265