

**Product Name :**  
Storage Vial

**Product Code :**  
CHEMISTRYUK-AS03LAB002R



**Description :**

Storage Vial

**Technical Specification :**

The Cap is provided with a built-in seal ring making the vial leakproof.

The base of the vial has longitudinal grooves which give it a non-twisting character when placed in a vial rack.

As the name suggests, these vials are commonly used for storage of biological material, human & animal cells at temperature as low as -190°C.

**Capacity:** 1.0ml, 2.0ml, 5.0ml



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265