

Product Name :
Storage Vial Internal Thread

Product Code :
CHEMISTRYSUK-AS03LAB004R



Description :

Storage Vial Internal Thread

Technical Specification :

The cap fitted with O-ring screw into the vial makes it a leakproof storage vial.

Another addition to the range of storage vials is this vial which has internal thread in it.

Capacity: 2.0ml

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific