

**Product Name :**  
SV10-SV5

**Product Code :**  
CHEMISTRYUK-AS03LAB005R



**Description :**

SV10-SV5

**Technical Specification :**

Available in 5ml and 10ml capacities, these vials have polypropylene cap with built in seal ring making it a leakproof vial.

These polypropylene vials add to the existing range of storage vials.

Just like others storage vials. These are also used for the storage of biological material, animal & human cells etc.

**Capacity:** SV-10ml, SV-5ml



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265