

Product Name :
Micro Tip Box

Product Code :
CHEMISTRYSUK-AS03LAB002S



Description :

Micro Tip Box

Technical Specification :

Micro Tip Box at most competitive prices. These Micro Tip Box available in different colour can easily fit the cap of Micro Tip Box for colour identification of different samples.

Capacity:.01 For 96 Micro Tips of 2-200 μ l

.02 For 100 Micro Tips of 200-1000 μ l

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

