

Product Name :
Separatory Funnel

Product Code :
CHEMISTRYSUK-AS03LAB003T



Description :

Separatory Funnel

Technical Specification :

These Separatory Funnel available in different colour can easily fit the cap of Separatory Funnel for colour identification of different samples.

Type Capacity (ml)

0.01 100

0.02 250

0.03 500



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265