

**Product Name :**  
Funnels Buchner

**Product Code :**  
CHEMISTRYSUK-AS03LAB004T



**Description :**

Funnels Buchner

**Technical Specification :**

These Funnels Buchner available in different colour can easily fit the cap of Funnels Buchner for colour identification of different samples.

Type Diameter (mm)

0.01 50

0.02 75

0.03 100

0.04 125

0.05 150



## Astra Scientific

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265