

**Product Name :**  
Funnels Buchner, Polypropylene

**Product Code :**  
CHEMISTRYSUK-AS03LAB005T



**Description :**

Funnels Buchner, Polypropylene

**Technical Specification :**

These Funnels Buchner, Polypropylene available in different colour can easily fit the cap of Funnels Buchner, Polypropylene for colour identification of different samples.

Type Capacity (mm)

0.01 70

0.02 110



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265