

Product Name :
Gas Generator,kipps Glass

Product Code :
CHEMISTRYSUK-AS03LAB004U



Description :

Gas Generator,kipps Glass

Technical Specification :

These Gas Generator,kipps available in different colour can easily fit the cap of Gas Generator,kipps for colour identification of different samples.

Type Capacity (ml)

0.01 250

0.02 500

0.03 1000



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265