

Product Name :
Spirit Lamp, Metal

Product Code :
CHEMISTRYSUK-AS03LAB003V



Description :

Spirit Lamp, Metal

Technical Specification :

These Spirit lamp, Metal available in different colour can easily fit the cap of Spirit lamp, Metal for colour identification of different samples.

Capacity (ml)

6001

0002



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265