

Product Name :
Test Tube Stand, 3 Tier

Product Code :
CHEMISTRYSUK-AS03LAB005W



Description :

Test Tube Stand, 3 Tier

Technical Specification :

These Test tube stand, 3 Tier available in different colour can easily fit the cap of Test tube stand, 3 Tier for colour identification of different samples.

Capacity

~~031m~~
x 31 Tubes

~~062m~~
x 31 Tubes

~~2003m~~
x 20 Tubes

~~254m~~
x 12 Tubes

~~255m~~
x 18 Tubes

~~326m~~
x 12 Tubes

~~257m~~
x 36 Tubes

~~038m~~
x 18 Tubes

~~059m~~
x 18 Tubes

~~081mm~~
x 12 Tubes

~~031m~~
x 62 Tubes

~~062m~~
x 62 Tubes

~~2013m~~
x 40 Tubes

Material

Polycarbonate

~~0314~~m
x 31 Tubes

~~0615~~m
x 31 Tubes

~~0016~~m
x 20 Tubes

~~0517~~m
x 12 Tubes

~~0318~~m
x 62 Tubes

~~0619~~m
x 62 Tubes

~~0020~~mm
x 40 Tubes

~~0521~~m
x 24 Tubes

~~0522~~m
x 18 Tubes

~~0523~~m
x 36 Tubes

~~0224~~m
x 12 Tubes

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific