

Product Name :
Conical Measures

Product Code :
CHEMISTRYUK-AS03LAB006F



Description :

Conical Measures

Technical Specification :

These Conical Measures are very good for measuring small quantities of reagents.

The Measures are clear, have good chemical resistance and have raised graduations for easy reading.

Moulded in Polypropylene, these autoclavable Measures have a hexagonal base which provides it excellent stability.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265