

Product Name :
Borosilicate Glass Burette

Product Code :
CHEMISTRYSUK-AS03LAB006H



Description :

Borosilicate Glass Burette

Technical Specification :

Graduated with pinchcock and jet.

Capacity:

25cm³ x 0.1cm³. Accuracy \pm 0.1cm³

50cm³ x 0.1cm³. Accuracy \pm 0.1cm³



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265