

**Product Name :**  
Burettes With Straight Bore PTFE (Teflon) Key Stopcock

**Product Code :**  
CHEMISTRYSUK-AS03LAB007H



**Description :**

Burettes With Straight Bore PTFE (Teflon) Key Stopcock

**Technical Specification :**

Borocilicate glass.

Accuracy as per class B.

**Capacity:-**

10ml

25ml

50ml



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265