

Product Name :
Retort Clamp Four Finger

Product Code :
CHEMISTRYSU-K-AS03LAB006M



Description :

Retort Clamp Four Finger

Technical Specification :

Die cast, strong for holding objects up to 100mm, corklined jaws provide secure contact, with planted brass tommy bar adjusting screw and aluminium anodized stem 150 x 8mm, length x dia. Without boss head.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265