

**Product Name :**  
Retort Clamp Three Finger

**Product Code :**  
CHEMISTRYSUK-AS03LAB007M



**Description :**

Retort Clamp Three Finger

**Technical Specification :**

Three Prong, die cast in alloy, with other parts of plated brass, clamp accepts articles up to 75mm . diameter, actuated by nickle plated brass pummels and twin screw. Jaws are covered in rubber. Fitted with 15cm. long mild steel or aluminium rod.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265