

**Product Name :**  
Retort Rings

**Product Code :**  
CHEMISTRYSUK-AS03LAB008M



**Description :**

Retort Rings

**Technical Specification :**

A Brass thumbscrew chrome plated to fit boss. Cast iron with boss, hammertone painted. Suitable for fixing on rods up to 13mm dia. The distance between center of ring and supporting rod is approx. 120mm for all sizes.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265