

**Product Name :**  
Test Tube Holder

**Product Code :**  
CHEMISTRYSUK-AS03LAB010M



**Description :**

Test Tube Holder

**Technical Specification :**

The jaws will accept tubes up to 38mm diameter and their alignment is maintained by a wire collar. Chrome plated spring steel wire holder with finger grips.

Overall length 130mm.

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com