

Product Name :
MCT Twin Rack

Product Code :
CHEMISTRYSUUK-AS03LAB007P



Description :

MCT Twin Rack

Technical Specification :

Capacity: Rack for 0.5ml and 1.5ml MCTs.

Provided this Twin Purpose Rack for both 0.5ml and 1.5ml micro centrifuge tubes.

Each side has 108 holes for 0.5 and 1.5ml MCTs, respectively.

Each face of the base has a unique alpha numeric index for easy identification of the tube.

The transparent top provides an excellent view as well as a protection cover for the tubes stored.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265