

Product Name :
Rack For Scintillation Vial

Product Code :
CHEMISTRYSUK-AS03LAB007R



Description :

Rack For Scintillation Vial

Technical Specification :

These rugged and stable racks can hold scintillation vials for convenient storage, handling & sample preparations.

These Racks for Scintillation Vials, moulded in polypropylene, are autoclavable & corrosion resistant.

There is an individual rack for both 8ml and 20ml vials.

Capacity:

Rack for 50 Vials of 20ml Cap.

Rack for 90 Vials of 8ml Cap.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265