

Product Name :
Cryo Vial Internal Thread

Product Code :
CHEMISTRYUK-AS03LAB009R



Description :

Cryo Vial Internal Thread

Technical Specification :

The cap fitted with an O-ring screw into the vials and makes it the leakproof Cryo Vial.

Adding to the existing range of Cryo Vials is this which has internal thread in it.

These vials are provided with white, prominent graduations and a long making area for sample identification.

Capacity: 2.0ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265